
Niut

USER'S MANUAL
ELECTRIC PUSSY

Product Overview

This product can be used every day, and the combination of physical weights and smart exercisers will help the pelvic floor muscles. It is recommended to keep the daily exercise time within **15-30 minutes** for better results.

Product Information

Material: Silicone/ABS

Charging time: 60 minutes


Function: 10 frequency vibration

Endurance: 70 minutes

Weight: 35.7g/69g/139.2g/104.3g

Waterproof grade: Whole body waterproof

Product size: 133*32mm / 164*32mm / 167*37mm

ON/OFF instructions


Turn on: Press and hold the button for three seconds until the indicator light turns on.

Switch modes: Short press to cycle through 10 vibration modes.

Turn off: Press and hold for three seconds to power off.

Switch to somatosensory mode:

Use the upper corrugated button.

FM mode: Click to activate frequency adjustment (FM).

Pause: Long press to pause.

Note: The remote control is not waterproof
The remote control does not have a power switch function.

Steps for Use

You can exercise in any of the following ways:

Physical Exercise

1. Wash the exerciser with warm water or antiseptic soap and dry it before the first use.
2. Apply lubricating oil to the exerciser and insert it 4-6 cm into the vagina by hand. Ensure the silicone drawstring remains outside the vagina for easy removal.
3. Engage in activities such as walking, climbing stairs, lifting, jumping, etc., while ensuring the exercise ball remains in place.

Smart Exercise

1. Use the **bottom button** on the remote control:
 - Long press to pause.
 - Single click to adjust frequency.
 - Use the upper **corrugated button** to switch to somatosensory mode.
2. After turning on the remote control, tap the **corrugated button** to enter somatosensory mode, which has 10 adjustable frequencies that vary based on step size.
3. Turn on the lower end key for remote control mode, which activates the exercise ball. Tap the remote control key to switch frequency or mode.

Note: The remote control does not control the on/off functionality of the sphere.

Charging Instructions

Sphere Charging:

- Use the pinhole charging port to connect the USB cable to a charger.
- For the first charge, allow **120 minutes**. Subsequent charges require **60 minutes**.

Remote Control Battery:

- The remote control uses a built-in button battery, which should be replaced every 6-12 months.
- If water accidentally enters the remote control, open the back, remove the battery, dry it thoroughly, and reassemble.

Storage Guidelines

- Clean and dry the product before storing it in a bag or box.
- Store in a dry, non-corrosive area away from direct sunlight and out of children's reach.

Safety Precautions

- This product is designed for adults aged **18 years and older**.
- Avoid exposure to fire, chemicals, flammable or explosive materials, and high temperatures.
- Do not share the product with others. If discomfort occurs during use, discontinue use immediately.

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